EXERCISE INSTRUCTION MANUAL

**THANK YOU
FOR PURCHASING XTEND, EXTENSION
EXERCISER WITH CONSTANT FORCE
TECHNOLOGY!**

The following exercises are not tailored to any specific sport or activity, but instead are designed to maintain and improve joint range of motion, dexterity, coordination, flexibility, strength, endurance, circulation and general health of your fingers, hand, wrist and/or forearm. XTend comes equipped hook and loop strips that can be used for extra-thick or swollen fingers as well as to use as a thumb adapter. For best results, exercise thumb independently from fingers, not while exercising other fingers at the same time.

Depending on the exercise and/or the resistance level, you may need to hold XTend with your non-exercising hand for additional stability.

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
INSTRUCTIONS:

Place hand thru the wide elastic band toward the backside of XTend so that the top of XTend lines up across the crease at the base of your fingers. Place each finger into a finger cup so the finger is snug, but comfortable.

NOTE: If fingers are too small, turn finger cup around so that the finger enters the narrow portion of the finger cup for better grip. If the small finger still slips, follow instructions for using one of the hook/loop strips to secure finger as noted.

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
1. FULL FINGER EXTENSION:

GOOD FOR:

Place hand thru the wide elastic band as noted on page 1 "Instructions". Use the non-exercising hand to hold and stabilize the XTend while slowly and fully extending all fingers at one time. Hold for 3-5 seconds. Relax and slowly return to starting position. Repeat slowly, as instructed by a healthcare professional.

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


2. ISOLATED FULL FINGER EXTENSION:

GOOD FOR:



Place hand thru the wide elastic band as noted on page 1 "Instructions". Extend one finger at a time, and hold for 3-5 seconds. Relax fingers and slowly return to starting position. Repeat slowly with each finger, as instructed by a healthcare professional.

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
Place hand thru the wide elastic band as noted on page 1 "Instructions". Then place fingers thru the thinner elastic band located on the top of XTend so that the joints being exercised are located directly in front of the thin elastic band as shown. Finally, place each finger into each finger cups so the finger is snug, but comfortable. While keeping the MP joints slightly flexed, extend only the PIP and DIP joints of the finger(s) you want to exercise and hold for 3-5 seconds. Relax fingers and slowly return to starting position. Repeat slowly with each finger, as instructed by a healthcare professional.

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
GOOD FOR:

Step 1:

Before placing fingers in finger cups, attempt exercise without finger cups in order to ensure proper technique. Place hand thru the wide elastic band as noted on page 1 "Instructions". While keeping the PIP and DIP joints flexed, extend only the MP joints (all fingers at the same time) and hold for 3-5 seconds. Relax and slowly return to starting position.

Step 2: Place each finger into a finger cup so the finger is snug, but comfortable. While keeping the PIP and DIP joints flexed, extend only the MP joints (all fingers at the same time) and hold for 3-5 seconds. Relax and slowly return to starting position. Repeat slowly with each finger, as instructed by a healthcare professional. For best results, stabilize XTend with the other hand.

NOTE:

If still having difficulty keeping finger cups on fingertips, use one of the hook/loop strips to secure the small finger. Another option would be to place the fingers between strips of thin elastic material on top of XTend, and while keeping the PIP and DIP joints flexed, extend only the MP joints (all fingers at the same time) and hold for 3-5 seconds. Relax and slowly return to starting position.

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6. ISOLATED FINGER ADDUCTION/ABDUCTION:

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Place hand thru the wide elastic band as noted on page 1 “Instructions”. Focusing on one finger at a time, extend each finger as far as possible, hold for 2 seconds, move the finger toward one side, hold for 2 seconds then move to other side and hold for 2 seconds and relax and slowly return to starting position. Repeat slowly for each finger, as instructed by a healthcare professional. For best results, stabilize XTend with the other hand.

NOTE:
For an additional challenge, place 1st finger in second cup, 2nd finger in 1st cup; 3rd finger in 4th cup, 4th finger in 3rd cup, so that retractable cords are “crisscrossed”. Now extend fingers and spread them apart, hold for 2 seconds and relax back to starting position.

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



Turn XTend upside down and place hand thru the wide elastic band toward the backside of XTend unit so that the bottom of the XTend product lines up with the crease at the base of your fingers. Place thumb in 1st finger cup. If thumb is too large for finger cup, see hook and loop adapter instructions. While keeping the base (CMC) of the thumb still, slowly extend the thumb IP and MP joints and hold for 3-5 seconds. Relax and slowly return to starting position. Repeat as instructed by a healthcare professional. For best results, stabilize XTend with the other hand.

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
8. THUMB EXTENSION/ ABDUCTION HURDLE:

GOOD FOR:

Turn XTend upside down and place hand thru the wide elastic band toward the backside of XTend unit so that the bottom of the XTend product lines up with the crease at the base of your fingers. Place thumb in 1st finger cup. If thumb is too large for finger cup, see hook and loop adapter instructions. Slowly extend the thumb IP and MP joints as noted in the Thumb Extension exercise, then move thumb out/around the base, away from the palm and fingers as if making a fist. Hold for 3-5 seconds. Relax and slowly return to starting position. Repeat as instructed by a healthcare professional. For best results, stabilize XTend with the other hand.

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9. THUMB TIP ISOLATED EXTENSION:

GOOD FOR:

Turn XTend upside down and place hand thru the wide elastic band toward the backside of XTend unit so that the bottom of the XTend product lines up with the crease at the base of your fingers. Place thumb through thinner band until the tip of the thumb (IP joint) is free to move, then, place thumb in 1st finger cup. If thumb is too large for finger cup, see hook and loop adapter instructions. Extend the only the thumb tip, hold for 3 seconds and relax. Repeat 5 times slowly as instructed by a healthcare professional.

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Caution:
LATEX
This product contains Natural Rubber Latex which may cause allergic reactions.

Designed in USA, made in China exclusively for Dynatomy.

D'Addario & Company, Inc,
595 Smith Street,
Farmingdale, NY 11735

Constant Force, Dynatomy and D'Addario are registered trademarks of D'Addario & Company, Inc., or its affiliates in the United States and/or other countries.

Copyright © 2018.
All Rights Reserved.
Patent Pending.

CEpartner 4U, Esdoornlaan 13
3951DB, Maarn, The Netherlands

Dynatomy™
(631) 815-3805
sales@dynatomyproducts.com

For more information, please visit dynatomyproducts.com